


# New Infinity 500 mL Bag Set with ENFit Connector



A new 500 mL enteral bag set for use with the Infinity Enteral Feeding Pump is now available:


- **What is different?** The new set comes with ONLY the ENFit connector, and is specifically for use with ENFit compatible feeding tubes.
- **What happens to the existing bag set?** The existing 500 mL bag set (INF0500-A) is still available for purchase, and will continue to include the ENFit connector system AND transitional stepped adapter.
- **What is the new order number?** INF0500-E
- **How are the bag sets differentiated?** Different colored labels are applied to the box. The existing bag set has a white label (below left), while the new ENFit-Only bag set has a teal label (below right).



## EXISTING INF0500-A

500 mL Enteral Pump Deliver Set with ENFit Connector System and Transitional Stepped Connector

**Unit of Measure:** Case  
**Quantity/Case:** 30 each



## NEW INF0500-E

500 mL Enteral Pump Delivery Set with ENFit Connector

**Unit of Measure:** Case  
**Quantity/Case:** 30 each



Speak live with a Moog Clinical Representative for pump questions and troubleshooting guidance 24 hours per day, 7 days per week.

Clinical & Customer Support **800.970.2337**  
Visit [infinityfeedingpump.com](http://infinityfeedingpump.com) for additional information

ENFit® is a registered trademark of GEDSA.  
Infinity® is a registered trademark of Zevox, Inc. Moog® is a registered trademark of Moog, Inc. © 2020 Moog, Inc. Moog Medical, 4314 Zevex Park Lane, Salt Lake City, UT 84123 USA. MD-81226 Rev. A

**MOOG®**  
Enhancing Healthcare.  
Enriching Lives.™